1. Produktbeschreibung

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1. Product Description

The Ventricular Drainage System consists of the following parts:

- Skull screw for attachment to the skull and introduction of the cannula.
- Cannula with a sharp stylet (to pierce the dura).
- Blunt stylet (for gentle advancement into the ventricle).
- Screwdriver (to insert the self-tapping skull screw)

The internal diameter of the cannula is 0.8 mm. Drainage is accomplished through 15 transverse bore holes with a diameter of 0.55 mm in an 18 mm area at the tip.

2. Intended Purpose

The system is used for ventriculostomy, for example, in cases of acute hydrocephalus, intracerebral haemorrhage into the ventricular system, subarachnoid bleeding or space-occupying tumours in the posterior cranial fossa with obstructive hydrocephalus.

2.1. Indications

The system is used for ventriculostomy, for example, in cases of acute hydrocephalus, intracerebral haemorrhage into the ventricular system, subarachnoid bleeding or space-occupying tumours in the posterior cranial fossa with obstructive hydrocephalus.

The drainage system can be installed in a few minutes and because of its stable fixation in the skull, it can remain in place for up to 30 days. Connection to a drainage system or a connecting tube for an irrigation cannula is accomplished by means of a Luer Lock connection. Due to this closed system, flushing and suction are possible without difficulty. Blocked cannulae can be exchanged without any difficulty.

2.2. Contraindications

When used together with anticoagulants, there is increased risk of brain haemorrhage. Prior to use, the user must determine whether the planned procedure in combination with anticoagulants could lead to unacceptable risks for the patient and if necessary, refrain from using this product together with anticoagulants.

2.3. Intended Users

The Duisburg Needle Ventricular Drainage System must only be utilized by qualified personnel, i.e., the products may only be used by persons who have undergone relevant professional training or instruction by an experienced surgeon.

2.4. Patient Population

The Duisburg Needle Ventricular Drainage System is intended to be used in patients requiring ventriculostomy.

2.5. Undesirable Side Effect

Analysis of the clinical literature shows that the following possible side-effects and risks may occur when these products are used:

- Reaccumulation (recurrence)
- Haemorrhage (subdural, epidural and intracerebral), especially in patients with coagulation disorders
- Cerebrovascular accident
- Epidural haematoma
- Hydrocephalus
- Pneumocephalus
- Seizures
- Infection/skin irritation (preventive: careful intraoperative haemostasis of the incision)
- Blocking of the screw (warning: haematoma persistence, possible recurrence with resulting increase in intracranial pressure)
- Screw breakage (if subjected to lateral force, e.g. if the patient falls)
- Mortality risk: 1% - 3%

2.6. Warnings

The following warnings are applicable:

Verify the intraventricular location of the needle tip by means of CT.

STERILE EO

The products are supplied in both sterile and non-sterile condition. The variants are clearly differentiated by their product label. Non-sterile products must be removed from their packaging prior to use and cleaned, disinfected and sterilised by the user according to the specified procedures.

Only parts from WILDMedical GmbH may be used; do not use articles from other manufacturers.

Reprocessing of disposable products may lead to reduced performance or loss of functionality. Reuse of disposable products may lead to exposure to viral, bacterial, fungal or prionic pathogens. There are no permissible cleaning and sterilisation methods or recommendations available for this product. These products must not be cleaned, disinfected or sterilised.

Products labelled as "STERILE" may be used only if the packaging has not been damaged and the specified expiration date has not expired. Observe the notes on sterilisation in the next section.

Non-sterile products must be removed from their packaging prior to cleaning, disinfection and sterilization, since the packaging is not suited for cleaning, disinfection and sterilization!

Products are MR unsafe.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

2.7. PRINCIPLES OF OPERATION


Ventricular Drainage System "Duisburg Needle" is a device implicated in ventriculostomy, that involves creating an artificial opening/hole, the stoma, within a cerebral ventricle for the diversion or drainage of CSF, generally contained within the ventricular system. This diversion or drainage becomes necessary when the normal flow of the CSF is obstructed, resulting in an increased ICP.

The standard target point for the needle is the anterior horn of the lateral ventricle at the virtual intersection between the ipsilateral corner of the eye and the external auditory canal. In case of any suspicion of displacement or narrowing of the ventricle, an examination such as a CT is recommended before inserting the needle. After raising the head of the patient's bed to about 30°, the skin at the operative site in the fronto-dorsal area at the level of the coronal suture is thoroughly shaved and washed immediately before the operative procedure. After thorough disinfection/antisepsis of the skin around the operative field (adhere to exposure times), the area around the operative field is draped with sterile waterproof drapes.

Further, 2 ml of a local anesthetic is injected at a point 20 mm parasagittal and 20 mm anterior to the coronal suture, followed by making a 5 mm incision using a number 11 scalpel blade (to prevent skin contusion when inserting the screw). After drilling a 1.5 to 2.0 mm hole, the self-tapping skull screw is inserted using the special screwdriver. The correct direction for drilling and the screw is the virtual intersection between the ipsilateral corner of the eye and the external auditory canal. Before inserting the cannula in the skull screw, it must be ensured that the lateral locking screw has been unscrewed sufficiently to allow the cannula to be inserted entirely as far as it will go (the intracranial length of the cannula is about 25 mm shorter than the nominal length). The cone of the cannula must sit snugly in the skull screw to ensure an optimal seal. The Luer lock attachment is then locked into the skull screw by tightening the lateral fixation screw to prevent the cannula from slipping out. The cannula is then inserted using the sharp mandrin. After penetrating the dura, the sharp mandrin is pulled back by about 2 mm and replaced with the blunt mandrin. There are markings at one-centimeter intervals on the cannula starting 4 cm from the tip, which assist with the orientation. Once the fluid flow has been initiated (it is possible to irrigate with 1 ml of saline solution and then carefully aspirate), a drainage system or an ICP sensor can be attached. The sterility should be maintained at all times, and the connection should be monitored. Antibiotic prophylaxis is not usually necessary; however, daily laboratory testing should be performed on the fluid and the screw seating. The loose screws should be hand tightened. If secure seating of the screw cannot be achieved, then the drain should be removed due to the risk of infection. The blocked cannulae can be exchanged without difficulty.

Cautions during operation:

- The skull screw or the cannula should not be bent.
- The screw should not be turned too far.
- The height of the system overflow should constantly be monitored.
- Since mechanical damage can lead to the breakage of the screw and the cannula in the skull, the restless patients should be monitored.

Functional Description: Insert the self-tapping skull screw using the special screwdriver. After screw fixation in the skull, the cannula is applied. This is guided through the cranial screw and (1-6) connected to a conventional drainage system.

2.8. Shelf Life (sterile products)

The shelf life of the device is 5 years. After 5 years of aging the packages stayed tight.

2.9. Basic-UDI-DI

The basic UDI-DI for the Ventricular Drainage System – Duisburg Needle is: 426235778MW0120XBD.

3. Intended Clinical Benefits

- Allows drainage of and access to cerebrospinal fluid.
- Allows for monitoring of intracranial pressure.

<https://ec.europa.eu/tools/eudamed>

"This is the Summary of Safety and Clinical Performance (SSCP) location after the launch of European Database on Medical Devices/ Eudamed".

4. Preparation and Sterilisation

Products are intended for single use only and must not be reprocessed!

STERILE EO

The Duisburg needle and its components are available in both a sterile and a non-sterile variant. Products labelled as "NON STERILE" may be cleaned, disinfected and sterilized according to the following instructions.

Before cleaning, disinfection and sterilisation, the products should be removed from their packaging, since the packaging is not suited for cleaning, disinfection and sterilisation.

The manufacturer has validated that the following instructions for the preparation of the products with regard to their cleaning, disinfection and sterilisation are SUITABLE. The user is responsible for ensuring that the cleaning, disinfection and sterilisation that are actually performed achieve the desired results at the preparation facility using its equipment, materials and personnel. To this end, validation and regular monitoring of the process are routinely required. Likewise, any deviation from the instructions provided should be thoroughly evaluated with respect to its efficacy and possible harmful consequences.

4.1. Alkaline automatic cleaning/disinfection

For automatic cleaning and disinfection, a washer-disinfector that complies with ISO 15883 must always be used.

Pre-cleaning:

Place instruments in cold water (tap water <40°C with drinking water quality) for at least five (5) minutes. If possible, disassemble the instruments and clean under cold water (tap water <40°C with drinking water quality) using a soft brush until no more residue is visible. Flush cavities, threads, and drill holes for at least ten (10) seconds using water (tap water <40°C with drinking water quality) from a nozzle attachment that is operated by quick switching on and off in a pulsed flow, and rebrush. This process is performed automatically in some machines. The process is used to remove deposits that may have penetrated this area.

Place instruments in an ultrasound bath for 10 minutes at 40°C with 0.5 percent alkaline cleaning agent and apply ultrasound.

Rinse instruments with cold water (tap water <40°C with drinking water quality).

Cleaning:

1. Place the instruments into a tray on the instrument rack and start the cleaning process.
2. Rinsing with cold water for one (1) minute
3. Draining
4. Rinsing with cold water for three (3) minutes
5. Draining
6. Washing at 55 C (+/- 2°C) with 0.5 percent alkaline cleaning agent for five (5) minutes
7. Draining
8. Neutralisation with hot tap water (>40°C) and neutralizer for three (3) minutes
9. Draining
10. Second rinse with deionised warm water (>40°C) for two (2) minutes
11. Draining

Thermal Disinfection:

Automatic thermal disinfection in the washer-disinfector taking into account the national requirements for the A0 value (see ISO 15883), e.g. A0 value 3000: > 90°C for five (5) minutes.

Drying

Automatic drying according to the automatic drying process of the washer-disinfector:

30 minutes at 90°C (+/- 2°C)

If necessary, subsequent manual drying with lint-free cloth and blowing out of lumens using sterile, oil-free compressed air.

Information on the validation of the preparation process:

The following test instructions, materials and machinery or equivalents were used for validation:

Cleaning agent:

- Neodisher Mediclean forte

Neutralising agent:

- Neodisher Z

Cleaning/disinfection device:

- Miele G 7835 CD (Des-Var-TD)

Instrument rack:

- Miele MIC rack E 450/1

4.2. Control and Inspection

After successful cleaning and disinfection procedures, the products should be examined for any damage, discolouration or deformation. Damaged, discoloured or deformed products should be discarded.

4.3. Packaging

For sterile packing, packaging suitable for sterilisation with moist heat (standard polyethylene/Tyvek pouch according to DIN EN ISO 11607 and EN 868) of sufficient size should be used. The packaging must not be under any tension.

4.4. Sterilisation

The products packaged in separate sterile packaging were sterilized in the fractionated pre-vacuum process according to DIN EN ISO 17665, taking into account applicable national requirements.

- Three pre-vacuum phases
- Heating to a sterilisation temperature of 134°C (273.2 F)
- Hold time 5 min.
- Drying time 20 min.

Sterilisation method: Flash sterilisation is not suitable for products with a lumen.

4.5. Storage

 For best storage results, store products in a dry, dark, and cool place. The rooms in which sterile products are stored must be easy to clean and should not be accessible to the general public.

4.6. After Use

With the exception of the screwdriver, the products are specified for single use. Used products must therefore be discarded. Products that are known to have been inserted in patients with CJD (Creutzfeldt-Jakob disease) or vCJD (variant Creutzfeldt-Jakob disease) or patients suspected of having these illnesses should be disposed of by incineration.

5. Symbols/Explanations

The symbols depicted on the product or product label have the following meanings:

	Symbol for "Caution"
	Symbol for "Non sterile"
	Symbol for "Sterilized using ethylene oxide"
	Symbol for "Do not re-use"
	Symbol for "Consult instructions for use"
	Symbol for "Date of manufacture"
	Symbol for "Use by date"
	Symbol for "Keep dry"
	Symbol for "Keep away from sunlight"
	Symbol for "Manufacturer"
	Symbol for "Lot number"
	Symbol for "Catalogue number"

	Symbol for "Do not use if package is damaged"
	CE marking with Notified Body identification number (DQS Medizinprodukte GmbH)
	Symbol for "CE Marking of Conformity"
	Symbol for "Medical Device"
	Symbol for "MR Unsafe"
	Symbol for "Single sterile barrier system"
	Symbol for "Single sterile barrier system inside protective packaging"

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Marcado CE con número de identificación del organismo notificado (DQS Medizinprodukte GmbH)
WILDmedical GmbH