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Instruction for use WM-TF03 - Disposable Bioptome

1. Product Description

Disposable bioptomes are suitable for multiple biopsies during a single procedure on the same patient. Disposable bioptomes are available in various lengths (510 mm to 1200 mm) and diameters (Ø1.6 mm, Ø1.8 mm, Ø2.2 mm). They are supplied sterile. **STERILE EO**

2. Intended Use

2.1. Indications

Disposable bioptomes are used exclusively for the removal of tissue samples from the heart (endoscopic myocardial biopsy) so that complications, such as rejection reactions after transplant, can be identified at an early stage.

2.2. Contraindications

Disposable bioptomes should not be used if aneurysms are present in the area of use.

2.3. Intended Clinical Benefits

- Facilitates multiple biopsies in a single procedure
- Biopsy sampling to aid in diagnosis of conditions such as unclear cardiomyopathy, myocarditis or unexplained ventricular tachyarrhythmias

<https://ec.europa.eu/tools/eudamed>

"This is the Summary of Safety and Clinical Performance (SSCP) location after the launch of European Database on Medical Devices/ Eudamed"

2.4. Target Population

The Bioptome device is used in patients requiring tissue biopsies from the heart.

2.5. Basis-UDI-DI

The basic UDI-DI for the Disposable Bioptomes is:

426235778MW011008S.

3. USE

The most important product features should be inspected by the user before use of the product, even though we have conducted a final inspection of the same kind ourselves. Dispose of damaged products.

THESE INSTRUCTIONS DO NOT REPLACE INSTRUCTION PROVIDED BY AN EXPERIENCED SURGEON! ONLY PERSONS WHO HAVE RECEIVED SPECIAL TRAINING OR INSTRUCTION MAY USE THE PRODUCTS.

To provide access, first position a sheath of appropriate size according to instructions. The heart can be accessed by inserting an appropriate guide catheter through the internal jugular vein, femoral vein, or subclavian vein. For this purpose, the surgeon has a choice of bioptomes available in different lengths and diameters. Provided the bioptome functions properly, it is introduced through the sheath into the right ventricle under ECG monitoring. The bioptome must be kept closed! In the event of persistent arrhythmias, the instrument should be removed. Advance the bioptome to the desired tissue sampling site under ultrasound or image converter monitoring—following radiation protection guidelines—and obtain the biopsy. Avoid taking multiple biopsies from the same site. After tissue sampling, carefully pull the bioptome back through the sheath, with the jaws shut. Then, open the jaws and place the biopsy in a petri dish filled with saline solution.

4. Information and Warnings

The shaft of the disposable bioptome can be bent gently by hand for accurate sampling. Bend the shaft at most up to 20 mm proximal to its working end to prevent the jaws from opening asymmetrically.

Check the instrument again once it has been performed as needed. Do not use the instrument if the mechanism cannot easily be operated.

Do not open or close disposable bioptomes past the stop. Forceful operation can destroy the mechanism.

Following the biopsy, the patient must be monitored for at least one hour, since possible latent injury might exist.

Products labelled as "STERILE" should only be used prior to the expiration date indicated and if the packaging is undamaged. Because of their design, they are not suitable for reprocessing or resterilisation!

Reprocessing of disposable products may lead to reduced performance or loss of functionality. Reuse of disposable products may lead to exposure to viral, bacterial, fungal or prionic pathogens. No permissible cleaning and sterilisation methods or recommendations for reprocessing to original specifications are available for this product. This product must not be cleaned, disinfected or sterilised.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

4.1. Possible Side-effects

Analysis of the clinical literature has revealed that the following possible side-effects and risks may result from the use of this product:

- blood vessel injuries,
- bruising/haematomas
- arrhythmias,
- possible ventricular tamponade,
- possible perforation of the free wall of the right ventricle,
- Introduction of bacteria into the bloodstream or heart muscle due to insufficient hygiene measure,
- Cardiac arrest
- Conduction abnormalities
- Possible right ventricular perforation with cardiac tamponade.
- Pneumothorax / Hemothorax
- Traumatic bicuspid valve injury due to repeated endomyocardial biopsy
- Pericardial effusion

4.2. Storage

For best storage results, store products in a dry, dark, and cool place. The rooms in which sterile products are stored must be easy to clean and should not be accessible to the general public.

5. Handling

All surgical instruments should always be handled with utmost care during transport and storage. This applies, in particular, to delicate tips and other sensitive areas.

6. Warranty

WILDmedical GmbH supplies to its customers only tested products that are free of defects.

All products are designed and manufactured to meet the highest quality standards. We will not accept liability for products that have been modified from the original, used for purposes other than their intended use, or treated or used improperly.

7. Symbols/Explanations

The symbols depicted on the product or product label have the following meanings:

	Symbol for "Caution"
	Symbol for "Sterilized using Ethylene Oxide"
	Symbol for "Do not Reuse"
	Symbol for "Consult Instructions for Use"
	Symbol for "Use By"
	Symbol for "Keep Dry"
	Symbol for "Keep Away From Sunlight"
	Symbol for "Manufacturer"
	Symbol for "Date of manufacture"
	Symbol for "Lot Number"
	Symbol for "Catalogue Number"
	Symbol for "Do Not Use if Package is Damaged"
	Symbol for "CE Mark - Symbol and Notified Body Identification"
	Symbol for "Medical Device"
	Symbol for "Single sterile barrier system"
	Symbol for "Single sterile barrier system inside protective packaging"

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Druckerei Hohl GmbH & Co. KG